

Product Datasheet

Human TRAM1 ELISA Kit (orb442725)

Description	Human Translocation Associated Membrane Protein 1 (TRAM1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q15629
Sensitivity	0.059ng/mL
Entrez	23471
Expiration Date	Please enquire.