

Product Datasheet

Porcine TREM1 ELISA Kit (orb442729)

Description	Porcine Triggering Receptor Expressed On Myeloid Cells 1 (TREM1) Sandwich ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q6TYI6
Sensitivity	11.7pg/mL
Entrez	396601
Expiration Date	Please enquire.