

Product Datasheet

Mouse TREM2 ELISA Kit (orb442731)

Description	Mouse Triggering Receptor Expressed On Myeloid Cells 2 (TREM2) Sandwich ELISA Kit
Reactivity	Mouse
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q99NH8
Sensitivity	23.8pg/mL
Entrez	83433
Expiration Date	Please enquire.