

Product Datasheet

Canine TRF ELISA Kit (orb442732)

Description	Canine Transferrin (TRF) Sandwich ELISA Kit
Reactivity	Canine
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Sensitivity	1.12ng/mL
Expiration Date	Please enquire.