

Product Datasheet

Equine TRF ELISA Kit (orb442733)

Description	Equine Transferrin (TRF) Sandwich ELISA Kit
Reactivity	Equine
Range	78.125-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P27425
Sensitivity	34ng/mL
Entrez	100146633
Expiration Date	Please enquire.