

Product Datasheet

Mouse TRF ELISA Kit (orb442735)

Description	Mouse Transferrin (TRF) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	7.812-500ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q921I1
Sensitivity	2.64ng/mL
Entrez	22041
Expiration Date	Please enquire.