

Product Datasheet

Rat TRH ELISA Kit (orb442740)

Description	Rat Thyrotropin Releasing Hormone (TRH) competitive inhibition ELISA Kit
Reactivity	Rat
Range	7.812-500pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P01150
Sensitivity	2.64pg/mL
Entrez	25569
Expiration Date	Please enquire.