

Product Datasheet

Mouse TRIP6 ELISA Kit (orb442741)

Description	Mouse Thyroid Hormone Receptor Interactor 6 (TRIP6) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9Z1Y4
Sensitivity	0.144ng/mL
Entrez	22051
Expiration Date	Please enquire.