

Product Datasheet

Human TRPM1 ELISA Kit (orb442744)

Description	Human Transient Receptor Potential Cation Channel Subfamily M, Member 1 (TRPM1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q7Z4N2
Sensitivity	0.051ng/mL
Entrez	4308
Expiration Date	Please enquire.