

Product Datasheet

Human TSG101 ELISA Kit (orb442758)

Description	Human Tumor Susceptibility Gene 101 (TSG101) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q99816
Sensitivity	0.105ng/mL
Entrez	7251
Expiration Date	Please enquire.