

Product Datasheet

Human TSPAN30 ELISA Kit (orb442765)

Description	Human Tetraspanin 30 (TSPAN30) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P08962
Sensitivity	0.058ng/mL
Entrez	967
Expiration Date	Please enquire.