

Product Datasheet

Mouse TSPAN30 ELISA Kit (orb442766)

Description	Mouse Tetraspanin 30 (TSPAN30) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma, tissue homogenates, cell lysates and other biological fluids.
Sample Types	plasma, tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P41731
Sensitivity	0.064ng/mL
Entrez	12512
Expiration Date	Please enquire.