

Product Datasheet

Mouse TTPa ELISA Kit (orb442771)

Description	Mouse Alpha Tocopherol Transfer Protein (TTPa) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q8BWP5
Sensitivity	0.055ng/mL
Entrez	50500
Expiration Date	Please enquire.