

Product Datasheet

Mouse TUBb1 ELISA Kit (orb442778)

Description	Mouse Tubulin Beta 1 (TUBb1) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	A2AQ07
Sensitivity	0.061ng/mL
Entrez	545486
Expiration Date	Please enquire.