

Product Datasheet

Human TWIST ELISA Kit (orb442782)

Description	Human Twist Transcription Factor (TWIST) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q15672
Sensitivity	0.055ng/mL
Entrez	7291
Expiration Date	Please enquire.