

Product Datasheet

Bovine TYR ELISA Kit (orb442789)

Description	Bovine Tyrosinase (TYR) Sandwich ELISA Kit
Reactivity	Bovine
Range	78.12-5000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q8MIU0
Sensitivity	25ng/mL
Entrez	280951
Expiration Date	Please enquire.