

Product Datasheet

Rat TYR ELISA Kit (orb442791)

Description	Rat Tyrosinase (TYR) Sandwich ELISA Kit
Reactivity	Rat
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	D4A9G4
Sensitivity	6.7pg/mL
Entrez	308800
Expiration Date	Please enquire.