

Product Datasheet

Human UBA52 ELISA Kit (orb442793)

Description	Human Ubiquitin A 52 Residue Ribosomal Protein Fusion Product 1 (UBA52) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P62987
Sensitivity	0.116ng/mL
Entrez	7311
Expiration Date	Please enquire.