

Product Datasheet

Human UBE2A ELISA Kit (orb442799)

Description	Human Ubiquitin Conjugating Enzyme E2A (UBE2A) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P49459
Sensitivity	0.081ng/mL
Entrez	7319
Expiration Date	Please enquire.