

Product Datasheet

VAV1 (phospho-Tyr174) antibody (orb4428)

Description	VAV1 (phospho-Tyr174) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human VAV1 around the phosphorylation site of Tyr174 (EI(p-Y)ED)
Target	phospho-VAV1 (Tyr174)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	98 kDa
Uniprot ID	P15498
Expiration Date	12 months from date of receipt.