

Product Datasheet

Human UCP2 ELISA Kit (orb442826)

Description	Human Uncoupling Protein 2, Mitochondrial (UCP2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P55851
Sensitivity	0.051ng/mL
Entrez	7351
Expiration Date	Please enquire.