

Product Datasheet

Rat UCP2 ELISA Kit (orb442827)

Description	Rat Uncoupling Protein 2, Mitochondrial (UCP2) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P56500
Sensitivity	0.055ng/mL
Entrez	54315
Expiration Date	Please enquire.