

Product Datasheet

Mouse UMOD ELISA Kit (orb442838)

Description	Mouse Uromodulin (UMOD) Sandwich ELISA Kit
Reactivity	Mouse
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, urine and other biological fluids.
Sample Types	serum, plasma, urine and other biological fluids.
Uniprot ID	Q91X17
Sensitivity	1.14ng/mL
Entrez	22242
Expiration Date	Please enquire.