

Product Datasheet

Porcine uPA ELISA Kit (orb442844)

Description	Porcine Plasminogen Activator, Urokinase (uPA) Sandwich ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P04185
Sensitivity	0.054ng/mL
Entrez	396985
Expiration Date	Please enquire.