

Product Datasheet

Mouse uPAR ELISA Kit (orb442848)

Description	Mouse Plasminogen Activator, Urokinase Receptor (uPAR) Sandwich ELISA Kit
Reactivity	Mouse
Range	125-8000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P35456
Sensitivity	47pg/mL
Entrez	18793
Expiration Date	Please enquire.