

Product Datasheet

Mouse UQCRC2 ELISA Kit (orb442853)

Description	Mouse Ubiquinol Cytochrome C Reductase Core Protein II (UQCRC2) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q9DB77
Sensitivity	0.063ng/mL
Entrez	67003
Expiration Date	Please enquire.