

Product Datasheet

Human USP14 ELISA Kit (orb442857)

Description	Human Ubiquitin Specific Peptidase 14 (USP14) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P54578
Sensitivity	0.108ng/mL
Entrez	9097
Expiration Date	Please enquire.