

Product Datasheet

Human USP6NL ELISA Kit (orb442858)

Description	Human USP6 N-Terminal Like Protein (USP6NL) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q92738
Sensitivity	0.058ng/mL
Entrez	9712
Expiration Date	Please enquire.