


Product Datasheet

Human Vaspin ELISA Kit (orb442869)

Description	Human Visceral Adipose Tissue Derived Serine Protease Inhibitor (Vaspin) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q8IW75
Sensitivity	0.054ng/mL
Entrez	145264
Expiration Date	Please enquire.



Standard Curve of Human Vaspin ELISA Kit...