

Product Datasheet

Porcine VCAM1 ELISA Kit (orb442877)

Description	Porcine Vascular Cell Adhesion Molecule 1 (VCAM1) Sandwich ELISA Kit
Reactivity	Porcine
Range	15.625-1000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q28939
Sensitivity	6.7ng/mL
Entrez	396925
Expiration Date	Please enquire.