

Product Datasheet

Rabbit VCAM1 ELISA Kit (orb442879)

Description	Rabbit Vascular Cell Adhesion Molecule 1 (VCAM1) Sandwich ELISA Kit
Reactivity	Rabbit
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	G1T1S9
Sensitivity	5.4pg/mL
Entrez	100008901
Expiration Date	Please enquire.