

Product Datasheet

Monkey VCAM1 ELISA Kit (orb442880)

Description	Monkey Vascular Cell Adhesion Molecule 1 (VCAM1) Sandwich ELISA Kit
Reactivity	Monkey
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q58T13
Sensitivity	0.36ng/mL
Entrez	712036
Expiration Date	Please enquire.