

Product Datasheet

Human VCP ELISA Kit (orb442886)

Description	Human Valosin Containing Protein (VCP) Sandwich ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P55072
Sensitivity	0.122ng/mL
Entrez	7415
Expiration Date	Please enquire.