

Product Datasheet

Mouse Visfatin ELISA Kit (orb442895)

Description	Mouse Visfatin (VF) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q99KQ4
Sensitivity	13.2pg/mL
Entrez	59027
Expiration Date	Please enquire.