

Product Datasheet

Rat vHL ELISA Kit (orb442898)

Description	Rat Von Hippel Lindau Tumor Suppressor (vHL) Sandwich ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	Q64259
Sensitivity	0.113ng/mL
Entrez	24874
Expiration Date	Please enquire.