

Product Datasheet

Human VIPR1 ELISA Kit (orb442907)

Description	Human Vasoactive Intestinal Peptide Receptor 1 (VIPR1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P32241
Sensitivity	0.045ng/mL
Entrez	7433
Expiration Date	Please enquire.