

Product Datasheet

Human Versican ELISA Kit (orb442920)

Description	Human Versican (VS) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q86W61
Sensitivity	27pg/mL
Entrez	1462
Expiration Date	Please enquire.