

Product Datasheet

Human VSNL1 ELISA Kit (orb442923)

Description	Human Visinin Like Protein 1 (VSNL1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cerebrospinal fluid and other biological fluids.
Sample Types	tissue homogenates, cerebrospinal fluid and other biological fluids.
Uniprot ID	P37235
Sensitivity	0.058ng/mL
Entrez	3241
Expiration Date	Please enquire.