

Product Datasheet

Human VTCN1 ELISA Kit (orb442924)

Description	Human V-Set Domain Containing T-Cell Activation Inhibitor 1 (VTCN1) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	Q7Z7D3
Sensitivity	0.054ng/mL
Entrez	79679
Expiration Date	Please enquire.