

Product Datasheet

Human vWA3A ELISA Kit (orb442930)

Description	Human Von Willebrand Factor A Domain Containing Protein 3A (vWA3A) Sandwich ELISA Kit
Reactivity	Human
Range	0.468-30ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	A6NCI4
Sensitivity	0.15ng/mL
Entrez	146177
Expiration Date	Please enquire.