

Product Datasheet

Porcine vWF ELISA Kit (orb442934)

Description	Porcine Von Willebrand Factor (vWF) competitive inhibition ELISA Kit
Reactivity	Porcine
Range	15.625-1000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma and other biological fluids.
Sample Types	plasma and other biological fluids.
Uniprot ID	Q28833
Sensitivity	4.51ng/mL
Entrez	399543
Expiration Date	Please enquire.