

Product Datasheet

Rat vWF ELISA Kit (orb442935)

Description	Rat Von Willebrand Factor (vWF) competitive inhibition ELISA Kit
Reactivity	Rat
Range	15.625-1000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	plasma.
Sample Types	plasma.
Uniprot ID	Q62935
Sensitivity	6.6ng/mL
Entrez	116669
Expiration Date	Please enquire.