

Product Datasheet

Human WNT11 ELISA Kit (orb442950)

Description	Human Wingless Type MMTV Integration Site Family, Member 11 (WNT11) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O96014
Sensitivity	0.058ng/mL
Entrez	7481
Expiration Date	Please enquire.