

Product Datasheet

Human WNT2 ELISA Kit (orb442952)

Description	Human Wingless Type MMTV Integration Site Family, Member 2 (WNT2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	P09544
Sensitivity	0.062ng/mL
Entrez	7472
Expiration Date	Please enquire.