

## Product Datasheet

### Mouse WNT3A ELISA Kit (orb442954)

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| <b>Description</b>       | Mouse Wingless Type MMTV Integration Site Family, Member 3A (WNT3A) Sandwich ELISA Kit   |
| <b>Reactivity</b>        | Mouse                                                                                    |
| <b>Range</b>             | 15.625-1000pg/mL                                                                         |
| <b>Form/Appearance</b>   | Traditional                                                                              |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only                                                                    |
| <b>Application notes</b> | tissue homogenates, cell lysates, cell culture supernates and other biological fluids.   |
| <b>Sample Types</b>      | tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| <b>Uniprot ID</b>        | <b>P27467</b>                                                                            |
| <b>Sensitivity</b>       | 6.7pg/mL                                                                                 |
| <b>Entrez</b>            | <b>22416</b>                                                                             |
| <b>Expiration Date</b>   | Please enquire.                                                                          |