

Product Datasheet

Rat WNT5A ELISA Kit (orb442963)

Description	Rat Wingless Type MMTV Integration Site Family, Member 5A (WNT5A) Sandwich ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	F2Z3U1
Sensitivity	0.055ng/mL
Entrez	64566
Expiration Date	Please enquire.