

Product Datasheet

Human WWP2 ELISA Kit (orb442969)

Description	Human WW Domain Containing E3 Ubiquitin Protein Ligase 2 (WWP2) Sandwich ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates, cell lysates and other biological fluids.
Sample Types	tissue homogenates, cell lysates and other biological fluids.
Uniprot ID	000308
Sensitivity	0.27ng/mL
Entrez	11060
Expiration Date	Please enquire.