

Product Datasheet

Human XCL2 ELISA Kit (orb442972)

Description	Human Chemokine C-Motif Ligand 2 (XCL2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9UBD3
Sensitivity	0.058ng/mL
Entrez	6846
Expiration Date	Please enquire.