

Product Datasheet

Bovine XDH ELISA Kit (orb442973)

Description	Bovine Xanthine Dehydrogenase (XDH) Sandwich ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P80457
Sensitivity	6.3pg/mL
Entrez	280960
Expiration Date	Please enquire.