

Product Datasheet

Mouse XRCC5 ELISA Kit (orb442984)

Description	Mouse X-Ray Repair Cross Complementing 5 (XRCC5) Sandwich ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	P27641
Sensitivity	0.072ng/mL
Entrez	22596
Expiration Date	Please enquire.